

MaternalCare - Doctor's Portal Test Documentation

Module: Alert System Test Cases

Test Case ID:	TC_DOC_ALERT_01		
Designed by:	Danny D. Dinglasa Jr.	Tested by:	Pheinz Reneil T. Magnun
Test Title:	Verify Real-time Alert Monitoring and Detail View		
Pre-conditions:	1. User is logged in as a Doctor. 2. There are active alerts in the system.		

#	Test Steps	Expected Result	Status	Recommendation
1	Observe the 'Alerts' tab in the bottom navigation bar.	A numbered badge appears if there are high/medium alerts. Red for high, Yellow for medium.	PASSED	Use intuitive icons and color-coding (e.g., Red for high risk, Green for normal) consistently.
2	Navigate to the Alerts tab.	Alerts are displayed in descending order of urgency.	PASSED	Implement skeleton loaders to improve perceived performance during data fetching.
3	Click 'View Details' on a High Priority alert.	User is taken to the Alert Details screen showing graphs and risk factors.	PASSED	Use intuitive icons and color-coding (e.g., Red for high risk, Green for normal) consistently.
4	Verify the 'Systolic BP Trend' line chart.	Chart shows previous BP readings to visualize trends.	PASSED	Ensure all UI elements meet WCAG 2.1 Level AA contrast standards for better accessibility.
5	Click 'Perform Re-assessment' from Alert Details.	User is taken to Assessment form with Patient Name and Age pre-filled.	PASSED	Provide clear, non-technical error messages to guide users during form submission errors.

MaternalCare - Doctor's Portal Test Documentation

Module: Dashboard Test Cases

Test Case ID:	TC_DOC_DASH_01		
Designed by:	Danny D. Dinglasa Jr.	Tested by:	Pheinz Reneil T. Magnun
Test Title:	Verify Doctor Dashboard Analytics and Navigation		
Pre-conditions:	1. User is logged in as a Doctor. 2. User is on the Dashboard screen.		

#	Test Steps	Expected Result	Status	Recommendation
1	Observe the summary cards (Total Patients, High Risk, etc.).	Cards display accurate counts from the database.	PASSED	Add descriptive tooltips to complex clinical charts and AI insight parameters.
2	Verify the Pie Chart for risk breakdown.	Pie chart displays visual distribution of Low, Medium, and High risk patients.	PASSED	Maintain consistent typography and spacing (padding/margins) throughout the application.
3	Find 'Recent Assessments' and click 'View all'.	User is redirected to the Patient Records screen.	PASSED	Use intuitive icons and color-coding (e.g., Red for high risk, Green for normal) consistently.
4	Click on a specific patient row in the 'Recent Assessments' list.	User is taken to the Patient Record Detail screen for that specific patient.	PASSED	Ensure the layout is fully responsive across mobile, tablet, and desktop viewports.
5	Click 'One-Click Share Access' and select 'Scan QR Code'.	Camera opens and successfully scans a patient's QR code, redirecting to their records.	PASSED	Implement skeleton loaders to improve perceived performance during data fetching.
6	Click 'One-Click Share Access' and enter a valid 6-digit code manually.	System verifies the code and grants access to the patient's record detail.	PASSED	Maintain consistent typography and spacing (padding/margins) throughout the application.

MaternalCare - Doctor's Portal Test Documentation

Module: Model Insights Test Cases

Test Case ID:	TC_DOC_MODEL_01		
Designed by:	Danny D. Dinglasa Jr.	Tested by:	Pheinz Reneil T. Magnun
Test Title:	Verify AI Model Information and Explainability Info		
Pre-conditions:	1. User is logged in as a Doctor. 2. User is on the Model tab.		

#	Test Steps	Expected Result	Status	Recommendation
1	View the 'Random Forest Classifier' info card.	Displays detailed description of the machine learning ensemble method.	PASSED	Use intuitive icons and color-coding (e.g., Red for high risk, Green for normal) consistently.
2	Scroll to view feature importance or dataset info.	Information about training data and model versioning is visible.	PASSED	Implement skeleton loaders to improve perceived performance during data fetching.

MaternalCare - Doctor's Portal Test Documentation

Module: Patient Assessment Test Cases

Test Case ID:	TC_DOC_ASSESS_01		
Designed by:	Danny D. Dinglasa Jr.	Tested by:	Pheinz Reneil T. Magnun
Test Title:	Verify AI-Powered Risk Assessment Flow		
Pre-conditions:	1. User is logged in as a Doctor. 2. User is on the New Assessment screen.		

#	Test Steps	Expected Result	Status	Recommendation
1	Select 'New Patient' and enter details.	Input is accepted, validated, and the system proceeds to the next stage of the workflow.	PASSED	Ensure the layout is fully responsive across mobile, tablet, and desktop viewports.
2	Enter vital signs that indicate high risk.	AI model processes the input data and generates a risk probability with clinical explanations.	PASSED	Ensure the layout is fully responsive across mobile, tablet, and desktop viewports.
3	Click 'Save Report'.	Information is successfully committed to the database and reflected in the patient record.	PASSED	Ensure all UI elements meet WCAG 2.1 Level AA contrast standards for better accessibility.
4	Perform a re-assessment for the same patient with normal vitals.	Backend state and frontend UI are synchronized, showing the most recent data updates.	PASSED	Add descriptive tooltips to complex clinical charts and AI insight parameters.

MaternalCare - Doctor's Portal Test Documentation

Module: Patient Records Test Cases

Test Case ID:	TC_DOC_RECORDS_01		
Designed by:	Danny D. Dinglasa Jr.	Tested by:	Pheinz Reneil T. Magnun
Test Title:	Verify Patient History and Record Management		
Pre-conditions:	1. User is logged in as a Doctor. 2. User is on the Patient Records screen.		

#	Test Steps	Expected Result	Status	Recommendation
1	Search for a patient by name.	List filters to show only matching patients.	PASSED	Implement skeleton loaders to improve perceived performance during data fetching.
2	Click on a patient to view details.	Displays all historical assessments with dates and risk levels.	PASSED	Implement skeleton loaders to improve perceived performance during data fetching.
3	Open an assessment record from the history list.	Modal opens showing full vital details and AI recommendations.	PASSED	Implement clear visual focus indicators for keyboard-only navigation users.

MaternalCare - Doctor's Portal Test Documentation

Module: Risk Assessment Test Cases

Test Case ID:	TC_DOC_ASM		
Designed by:	Danny D. Dinglasa Jr.	Tested by:	Pheinz Reneil T. Magnun
Test Title:	Manual Verification: Risk Assessment		
Pre-conditions:	In Assessment form.		

#	Test Steps	Expected Result	Status	Recommendation
1	Switch temperature unit between °F and °C in assessment form.	UI updates placeholder and validation helpers to match selected unit.	PASSED	Provide clear, non-technical error messages to guide users during form submission errors.
2	Complete all 4 steps of the risk assessment form.	System calculates probability and generates clinical insights/recommendations.	PASSED	Use intuitive icons and color-coding (e.g., Red for high risk, Green for normal) consistently.
3	Submit assessment with Celsius input.	Fahrenheit conversion performed correctly before ML processing.	PASSED	Implement skeleton loaders to improve perceived performance during data fetching.

General System Recommendations

- Conduct quarterly penetration testing to identify and patch new vulnerabilities.
- Ensure all sensitive data is encrypted at rest using industry-standard AES-256.
- Implement multi-factor authentication (MFA) for administrative and medical staff accounts.
- Maintain a regular schedule for rotating API keys and database credentials every 90 days.
- Establish real-time monitoring and alerting for suspicious activity logs.
- Perform regular integrity checks on data backups and conduct disaster recovery drills.
- Ensure all third-party dependencies are audited and updated to their latest stable versions regularly.
- Implement a strict 'Least Privilege' access control model across all system modules.

MaternalCare - Doctor's Portal Test Documentation

Final Testing Summary

Test Case ID	Title	Remarks	Recommendation
TC_DOC_ALERT_01	Verify Real-time Alert Monitoring and Detail View	PASSED	Use intuitive icons and color-coding (e.g., Red for high risk, Green for normal) consistently.
TC_DOC_DASH_01	Verify Doctor Dashboard Analytics and Navigation	PASSED	Add descriptive tooltips to complex clinical charts and AI insight parameters.
TC_DOC_MODEL_01	Verify AI Model Information and Explainability Info	PASSED	Use intuitive icons and color-coding (e.g., Red for high risk, Green for normal) consistently.
TC_DOC_ASSESS_01	Verify AI-Powered Risk Assessment Flow	PASSED	Ensure the layout is fully responsive across mobile, tablet, and desktop viewports.
TC_DOC_RECORDS_01	Verify Patient History and Record Management	PASSED	Implement skeleton loaders to improve perceived performance during data fetching.
TC_DOC_ASM	Manual Verification: Risk Assessment	PASSED	Provide clear, non-technical error messages to guide users during form submission errors.

Conclusion: All functional test cases listed above have achieved a status of 'PASSED'. The verified UI workflows and AI-driven data retrieval demonstrate full compliance with the operational requirements for the Doctor's Portal.